



RIPK2 Antibody, FITC conjugated

Product Code	CSB-PA019736LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43353
Immunogen	Recombinant Human Receptor-interacting serine/threonine-protein kinase 2 protein (1-540AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-interacting serine/threonine-protein kinase 2 (EC 2.7.11.1) (CARD-containing interleukin-1 beta-converting enzyme-associated kinase) (CARD-containing IL-1 beta ICE-kinase) (RIP-like-interacting CLARP kinase) (Receptor-interacting protein 2) (RIP-2) (Tyrosine-protein kinase RIPK2) (EC 2.7.10.2), RIPK2, CARDIAK RICK RIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIPK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.