



SCP2 Antibody, Biotin conjugated

Product Code	CSB-PA020856DD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22307
Immunogen	Recombinant Human Non-specific lipid-transfer protein (1-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-specific lipid-transfer protein (NSL-TP) (EC 2.3.1.176) (Propanoyl-CoA C-acyltransferase) (SCP-chi) (SCPX) (Sterol carrier protein 2) (SCP-2) (Sterol carrier protein X) (SCP-X), SCP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SCP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.