



GRP Antibody, FITC conjugated

Product Code	CSB-PA009940EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07492
Immunogen	Recombinant Human Gastrin-releasing peptide protein (31-98AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrin-releasing peptide (GRP) [Cleaved into: Neuromedin-C (GRP-10)], GRP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRP