





PGLYRP2 Antibody

Product Code	CSB-PA846684LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PD5
Immunogen	Recombinant Human N-acetylmuramoyl-L-alanine amidase protein (22-576AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylmuramoyl-L-alanine amidase (EC 3.5.1.28) (Peptidoglycan recognition protein 2) (Peptidoglycan recognition protein long) (PGRP-L), PGLYRP2, PGLYRPL PGRPL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PGLYRP2

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA846684LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA846684LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.