



ACTA1 Antibody, FITC conjugated

Product Code	CSB-PA001205LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68133
Immunogen	Recombinant Human Actin, alpha skeletal muscle protein (3-377AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, alpha skeletal muscle (Alpha-actin-1) [Cleaved into: Actin, alpha skeletal muscle, intermediate form], ACTA1, ACTA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.