



aglB Antibody, Biotin conjugated

Product Code	CSB-PA381525LD01AVE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y865
Immunogen	Recombinant Aspergillus niger Probable alpha-galactosidase B protein (21-443AA)
Raised In	Rabbit
Species Reactivity	Aspergillus niger
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable alpha-galactosidase B (EC 3.2.1.22) (Melibiase B), aglB
Immunogen Species	Aspergillus niger
Research Area	Others
Target Names	aglB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.