



pntB Antibody, Biotin conjugated

Product Code	CSB-PA364617LD01EVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AB67
Immunogen	Recombinant Escherichia coli NAD(P) transhydrogenase subunit beta protein (330-462AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD(P) transhydrogenase subunit beta (EC 1.6.1.2) (Nicotinamide nucleotide transhydrogenase subunit beta) (Pyridine nucleotide transhydrogenase subunit beta), pntB
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	pntB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.