




MMP12 Antibody

Product Code	CSB-PA12379A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39900
Immunogen	Recombinant Human Macrophage metalloelastase protein (109-470AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage metalloelastase (MME) (EC 3.4.24.65) (Macrophage elastase) (ME) (hME) (Matrix metalloproteinase-12) (MMP-12), MMP12, HME
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP12


Image



Western blot
All lanes: MMP12 antibody at 2µg/ml
Lane 1: Hela whole cell lysate
Lane 2: NIH/3T3 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 55 kDa
Observed band size: 55 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA12379A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA12379A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.