



NAPSA Antibody, Biotin conjugated

Product Code	CSB-PA11529D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96009
Immunogen	Recombinant Human Napsin-A protein (64-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Napsin-A (EC 3.4.23.-) (Aspartyl protease 4) (ASP4) (Asp 4) (Napsin-1) (TA01/TA02), NAPSA, NAP1 NAPA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NAPSA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.