



QPCT Antibody, FITC conjugated

Product Code	CSB-PA019135LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16769
Immunogen	Recombinant Human Glutaminyl-peptide cyclotransferase protein (29-361AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaminyl-peptide cyclotransferase (EC 2.3.2.5) (Glutaminyl cyclase) (QC) (sQC) (Glutaminyl-tRNA cyclotransferase) (Glutamyl cyclase) (EC), QPCT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	QPCT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.