



RETN Antibody, HRP conjugated

Product Code	CSB-PA019573LB01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q762I5
Immunogen	Recombinant Bovine Resistin protein (19-109AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Resistin, RETN, RSTN
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	RETN