



TSPAN7 Antibody, Biotin conjugated

Product Code	CSB-PA025165LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41732
Immunogen	Recombinant Human Tetraspanin-7 protein (113-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-7 (Tspan-7) (Cell surface glycoprotein A15) (Membrane component chromosome X surface marker 1) (T-cell acute lymphoblastic leukemia-associated antigen 1) (TALLA-1) (Transmembrane 4 superfamily member 2) (CD antigen CD231), TSPAN7, A15 DXS1692E MXS1 TM4SF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSPAN7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.