



TSPAN7 Antibody, Biotin conjugated

Product Code	CSB-PA025165LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41732
Immunogen	Recombinant Human Tetraspanin-7 protein (113-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-7 (Tspan-7) (Cell surface glycoprotein A15) (Membrane component chromosome X surface marker 1) (T-cell acute lymphoblastic leukemia-associated antigen 1) (TALLA-1) (Transmembrane 4 superfamily member 2) (CD antigen CD231), TSPAN7, A15 DXS1692E MXS1 TM4SF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSPAN7