



DLK1 Antibody, FITC conjugated

Product Code	CSB-PA006945LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80370
Immunogen	Recombinant Human Protein delta homolog 1 protein (24-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	May have a role in neuroendocrine differentiation.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein delta homolog 1 (DLK-1) (pG2) [Cleaved into: Fetal antigen 1 (FA1)], DLK1, DLK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DLK1