



P Antibody, Biotin conjugated

Product Code	CSB-PA303154LD01RAH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69480
Immunogen	Recombinant Rabies virus Phosphoprotein protein (1-297AA)
Raised In	Rabbit
Species Reactivity	Rabies virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoprotein (Protein P) (Protein M1), P
Immunogen Species	Rabies virus
Research Area	Others
Target Names	P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.