





Hras Antibody

Product Code	CSB-PA15539A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20171
Immunogen	Recombinant Rat GTPase HRas protein (2-186AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse, Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase HRas (H-Ras-1) (Transforming protein p21) (c-H-ras) (p21ras) [Cleaved into: GTPase HRas, N-terminally processed], Hras, Hras1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Hras

Image



Western blot
 All lanes: Hras antibody at 2µg/ml
 Lane 1: 293T whole cell lysate
 Lane 2: Mouse brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 22 kDa
 Observed band size: 22 kDa



Immunohistochemistry of paraffin-embedded human lymphoid tissue using CSB-PA15539A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.