



Nid1 Antibody, Biotin conjugated

Product Code	CSB-PA015802LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10493
Immunogen	Recombinant Mouse Nidogen-1 protein (428-665AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nidogen-1 (NID-1) (Entactin), Nid1, Ent
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Nid1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.