




Nid1 Antibody

Product Code	CSB-PA015802LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10493
Immunogen	Recombinant Mouse Nidogen-1 protein (428-665AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nidogen-1 (NID-1) (Entactin), Nid1, Ent
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Nid1

Image



Western blot
 All lanes: Nid1 antibody at 6µg/ml + Rat lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 137 kDa
 Observed band size: 150 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.