




Nox4 Antibody

Product Code	CSB-PA015961LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q924V1
Immunogen	Recombinant Rat NADPH oxidase 4 protein (210-424AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH oxidase 4 (EC 1.6.3.-) (Kidney oxidase-1) (KOX-1) (Kidney superoxide-producing NADPH oxidase), Nox4, Kox
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Nox4

Image



Western Blot

Positive WB detected in: U251 whole cell lysate,

SH-SY5Y whole cell lysate

All lanes: Nox4 antibody at 3.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67 kDa

Observed band size: 67 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.