



cry3Aa Antibody, Biotin conjugated

Product Code	CSB-PA358337LD01BDG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A379
Immunogen	Recombinant Bacillus thuringiensis subsp. Tenebrionis Pesticidal crystal protein cry3Aa protein (58-644AA)
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. Tenebrionis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pesticidal crystal protein Cry3Aa (73 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIIIA(a)), cry3Aa, bt13 cry3A cryC cryIIa cryIIIA(a)
Immunogen Species	Bacillus thuringiensis subsp. tenebrionis
Research Area	Others
Target Names	cry3Aa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.