




Mup11 Antibody

Product Code	CSB-PA361419LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04938
Immunogen	Recombinant Mouse Major urinary protein 11 protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major urinary protein 11, Mup11, Mup9
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Mup11

Image



Western blot
 All lanes: Mup11 antibody at 2µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 21 kDa
 Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.