



sspB Antibody, Biotin conjugated

Product Code	CSB-PA859203ED01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99V46
Immunogen	Recombinant Staphylococcus aureus Staphopain B protein (220-393AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Staphopain B (EC 3.4.22.-) (Staphylococcal cysteine proteinase B) (Staphylopain B), sspB
Immunogen Species	Staphylococcus aureus (strain Mu50 / ATCC 700699)
Research Area	Others
Target Names	sspB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.