



# Idh1 Antibody

<b>Product Code</b>	CSB-PA348716EA01LMS
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P56512
<b>Immunogen</b>	Recombinant Lactobacillus plantarum L-lactate dehydrogenase 1 protein (1-320AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Lactobacillus plantarum
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	L-lactate dehydrogenase 1 (L-LDH 1) (EC 1.1.1.27), Idh1, I-IdhL Idh IdhL IdhL1
<b>Immunogen Species</b>	Lactobacillus plantarum (strain ATCC BAA-793 / NCIMB 8826 / WCFS1)
<b>Research Area</b>	Others
<b>Target Names</b>	Idh1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.