



Idh1 Antibody, HRP conjugated

Product Code	CSB-PA348716EB01LMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56512
Immunogen	Recombinant Lactobacillus plantarum L-lactate dehydrogenase 1 protein (1-320AA)
Raised In	Rabbit
Species Reactivity	Lactobacillus plantarum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-lactate dehydrogenase 1 (L-LDH 1) (EC 1.1.1.27), Idh1, I-IdhL Idh IdhL IdhL1
Immunogen Species	Lactobacillus plantarum (strain ATCC BAA-793 / NCIMB 8826 / WCFS1)
Research Area	Others
Target Names	Idh1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.