



ABI5 Antibody

Product Code	CSB-PA879765LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SJN0
Immunogen	Recombinant Arabidopsis thaliana Protein ABSCISIC ACID-INSENSITIVE 5 protein (1-442AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ABSCISIC ACID-INSENSITIVE 5 (Dc3 promoter-binding factor 1) (AtDPBF1) (Protein GROWTH-INSENSITIVITY TO ABA 1) (bZIP transcription factor 39) (AtbZIP39), ABI5, BZIP39 DPBF1 GIA1 NEM1
Immunogen Species	Arabidopsis thaliana
Research Area	Cell Biology
Target Names	ABI5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.