



TMEM25 Antibody, FITC conjugated

Product Code	CSB-PA023821LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YD3
Immunogen	Recombinant Human Transmembrane protein 25 protein (27-322AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 25, TMEM25
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.