



FSHB Antibody, Biotin conjugated

Product Code	CSB-PA12335D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01225
Immunogen	Recombinant Human Follitropin subunit beta protein (19-129AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Follitropin subunit beta (Follicle-stimulating hormone beta subunit) (FSH-B) (FSH-beta) (Follitropin beta chain), FSHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FSHB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.