




GEMIN7 Antibody

Product Code	CSB-PA862039EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H840
Immunogen	Recombinant Human Gem-associated protein 7 protein (1-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gem-associated protein 7 (Gemin-7) (SIP3), GEMIN7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GEMIN7

Image



Western blot

All lanes: GEMIN7 antibody at 2µg/ml

Lane 1: HeLa whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 15 kDa

Observed band size: 15 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.