



ILKAP Antibody, HRP conjugated

Product Code	CSB-PA011686HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0C8
Immunogen	Recombinant Human Integrin-linked kinase-associated serine/threonine phosphatase 2C protein (1-392AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin-linked kinase-associated serine/threonine phosphatase 2C (ILKAP) (EC 3.1.3.16), ILKAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ILKAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.