




ILKAP Antibody

Product Code	CSB-PA011686HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0C8
Immunogen	Recombinant Human Integrin-linked kinase-associated serine/threonine phosphatase 2C protein (1-392AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin-linked kinase-associated serine/threonine phosphatase 2C (ILKAP) (EC 3.1.3.16), ILKAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ILKAP

Image



Western blot

All lanes: ILKAP antibody at 2µg/ml

Lane 1: Rat kidney tissue

Lane 2: U251 whole cell lysate

Lane 3: Mouse stomach tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43 kDa

Observed band size: 43 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.