




HK1 Antibody

| | |
|----------------------------|-----------------------------------------------------------------------------------|
| Product Code | CSB-PA010468HA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19367 |
| Immunogen | Recombinant Human Hexokinase-1 protein (13-475AA) |
| Raised In | Rabbit |
| Species Reactivity | Human, Rat |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:1000-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Hexokinase-1 (EC 2.7.1.1) (Brain form hexokinase) (Hexokinase type I) (HK I), HK1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | HK1 |

Image



Western blot

All lanes: HK1 antibody at 2µg/ml + Rat brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 103, 102 kDa

Observed band size: 103 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.