



HK1 Antibody, HRP conjugated

Product Code	CSB-PA010468HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19367
Immunogen	Recombinant Human Hexokinase-1 protein (13-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hexokinase-1 (EC 2.7.1.1) (Brain form hexokinase) (Hexokinase type I) (HK I), HK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.