



ZRANB2 Antibody, HRP conjugated

Product Code	CSB-PA027127HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95218
Immunogen	Recombinant Human Zinc finger Ran-binding domain-containing protein 2 protein (1-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger Ran-binding domain-containing protein 2 (Zinc finger protein 265) (Zinc finger, splicing), ZRANB2, ZIS ZNF265
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZRANB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.