



# HNRNPA2B1 Antibody, Biotin conjugated

Product Code	CSB-PA010602DD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22626
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoproteins A2/B1 protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A2/B1), HNRNPA2B1, HNRPA2B1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HNRNPA2B1