



M Antibody, Biotin conjugated

Product Code	CSB-PA322192LD01RAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15200
Immunogen	Recombinant Rabies virus Matrix protein (1-202AA)
Raised In	Rabbit
Species Reactivity	Rabies virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein (Phosphoprotein M2), M
Immunogen Species	Rabies virus (strain PM1503/AVO1) (RABV)
Research Area	Others
Target Names	M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.