



# SLC27A2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA021534LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O14975
<b>Immunogen</b>	Recombinant Human Very long-chain acyl-CoA synthetase protein (283-620AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Very long-chain acyl-CoA synthetase (VLACS) (VLCS) (Fatty acid transport protein 2) (FATP-2) (Fatty-acid-coenzyme A ligase, very long-chain 1) (Long-chain-fatty-acid--CoA ligase) (Solute carrier family 27 member 2) (THCA-CoA ligase) (Very long-chain-fatty-acid-CoA ligase), SLC27A2, ACSVL1, FACVL1, FATP2, VLACS
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	SLC27A2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.