



ioIG Antibody, Biotin conjugated

Product Code	CSB-PA329675HD01BRG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26935
Immunogen	Recombinant Bacillus subtilis Inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase protein (1-344AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase (EC 1.1.1.18) (EC 1.1.1.369) (Myo-inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase) (MI 2-dehydrogenase/DCI 3-dehydrogenase), ioIG, idh
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	ioIG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.