



# pyrF Antibody, Biotin conjugated

Product Code	CSB-PA14337D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56155
Immunogen	Recombinant Helicobacter pylori Orotidine 5'-phosphate decarboxylase protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Orotidine 5'-phosphate decarboxylase (EC 4.1.1.23) (OMP decarboxylase) (OMPDCase) (OMPdecase), pyrF
Immunogen Species	Helicobacter pylori
Research Area	Others
Target Names	pyrF