



pyrF Antibody, Biotin conjugated

Product Code	CSB-PA14337D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56155
Immunogen	Recombinant Helicobacter pylori Orotidine 5'-phosphate decarboxylase protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Orotidine 5'-phosphate decarboxylase (EC 4.1.1.23) (OMP decarboxylase) (OMPDCase) (OMPdecase), pyrF
Immunogen Species	Helicobacter pylori
Research Area	Others
Target Names	pyrF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.