



cpg2 Antibody, HRP conjugated

Product Code	CSB-PA361895DB01PW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06621
Immunogen	Recombinant Pseudomonas sp.Carboxypeptidase G2 protein (23-415AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase G2 (CPDG2) (EC 3.4.17.11) (Folate hydrolase G2) (Glutamate carboxypeptidase) (Pteroylmonoglutamic acid hydrolase G2) (Glucarpidase), cpg2
Immunogen Species	Pseudomonas sp.
Research Area	Others
Target Names	cpg2