



DPP6 Antibody, Biotin conjugated

Product Code	CSB-PA007140YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42658
Immunogen	Recombinant Human Dipeptidyl aminopeptidase-like protein 6 protein (117-859AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl aminopeptidase-like protein 6 (DPPX) (Dipeptidyl aminopeptidase-related protein) (Dipeptidyl peptidase 6) (Dipeptidyl peptidase IV-like protein) (Dipeptidyl peptidase VI) (DPP VI), DPP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPP6