



DPP6 Antibody, Biotin conjugated

Product Code	CSB-PA007140YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42658
Immunogen	Recombinant Human Dipeptidyl aminopeptidase-like protein 6 protein (117-859AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl aminopeptidase-like protein 6 (DPPX) (Dipeptidyl aminopeptidase-related protein) (Dipeptidyl peptidase 6) (Dipeptidyl peptidase IV-like protein) (Dipeptidyl peptidase VI) (DPP VI), DPP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.