





ACTB Antibody

Product Code	CSB-PA001207NJ01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Peptide sequence from Human Actin, cytoplasmic 1 protein (102-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 1 (Beta-actin) [Cleaved into: Actin, cytoplasmic 1, N-terminally processed], ACTB
Immunogen Species	Homo sapiens (Human)
Research Area	Isotype/Loading Controls
Target Names	ACTB


Image



Immunohistochemistry analysis of human small intestine tissue using CSB-PA001207NJ01HU at dilution of 1:100



Immunohistochemistry analysis of human tonsil tissue using CSB-PA001207NJ01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA001207NJ01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.