





ACTB Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA001207NJ01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P60709 |
| Immunogen | Peptide sequence from Human Actin, cytoplasmic 1 protein (102-375AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Actin, cytoplasmic 1 (Beta-actin) [Cleaved into: Actin, cytoplasmic 1, N-terminally processed], ACTB |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Isotype/Loading Controls |
| Target Names | ACTB |


Image



Immunohistochemistry analysis of human small intestine tissue using CSB-PA001207NJ01HU at dilution of 1:100



Immunohistochemistry analysis of human tonsil tissue using CSB-PA001207NJ01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA001207NJ01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)