



cry1Ac Antibody, HRP conjugated

Product Code	CSB-PA356448EB01BDC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05068
Immunogen	Recombinant Bacillus thuringiensis subsp. kurstaki Pesticidal crystal protein Cry1Ac protein (972-1178AA)
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. kurstaki
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pesticidal crystal protein Cry1Ac (133 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIA(c)), cry1Ac, cry1A(c) cry218 cryIA(c)
Immunogen Species	Bacillus thuringiensis subsp. kurstaki
Research Area	Others
Target Names	cry1Ac
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.