




cry1Ac Antibody

Product Code	CSB-PA356448EA01BDC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05068
Immunogen	Recombinant Bacillus thuringiensis subsp. kurstaki Pesticidal crystal protein Cry1Ac protein (972-1178AA)
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. kurstaki
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus thuringiensis subsp. kurstaki
Research Area	Others
Target Names	cry1Ac

Image



Western Blot

Positive WB detected in: Recombinant cry1Ac protein at 20ng, 10ng, 5ng

All lanes: cry1Ac antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51 kDa

Observed band size: 51 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.