



FABP1 Antibody, Biotin conjugated

Product Code	CSB-PA007940YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07148
Immunogen	Recombinant Human Fatty acid-binding protein, liver protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein, liver (Fatty acid-binding protein 1) (Liver-type fatty acid-binding protein) (L-FABP), FABP1, FABPL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FABP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.