





GTF2H5 Antibody

Product Code	CSB-PA010014EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZYL4
Immunogen	Recombinant Human General transcription factor IIH subunit 5 protein (1-71AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	General transcription factor IIH subunit 5 (General transcription factor IIH polypeptide 5) (TFB5 ortholog) (TFIIH basal transcription factor complex TTD-A subunit), GTF2H5, C6orf175 TTDA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GTF2H5


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA010014EA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA010014EA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA010014EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.