



# luxS Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA541717HC01ENV   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P45578  |
| <b>Immunogen</b>           | Recombinant Escherichia coli S-ribosylhomocysteine lyase protein (2-171AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Escherichia coli  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | S-ribosylhomocysteine lyase (EC 4.4.1.21) (AI-2 synthesis protein)<br>(Autoinducer-2 production protein LuxS), luxS, ygaG |
| <b>Immunogen Species</b>   | Escherichia coli (strain K12)   |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | luxS  |