



luxS Antibody, HRP conjugated

Product Code	CSB-PA541717HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45578
Immunogen	Recombinant Escherichia coli S-ribosylhomocysteine lyase protein (2-171AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-ribosylhomocysteine lyase (EC 4.4.1.21) (AI-2 synthesis protein) (Autoinducer-2 production protein LuxS), luxS, ygaG
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	luxS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.